



ARRAY

7

NEUROLOGICAL AUTOIMMUNE REACTIVITY SCREEN™



2 mL Serum

Antigens Tested (IgG + IgA Combined)

- Myelin Basic Protein
- Asialoganglioside
- Alpha + Beta Tubulin
- Cerebellar
- Synapsin

CLINICAL USE:

- Assist in the early detection of neuro-autoimmunity.
- Evaluate severity of the autoimmune response.
- Assist in monitoring the effectiveness of related treatment protocols.

RECOMMENDED FOR PATIENTS WHO:

- Have idiopathic neurological conditions such as Cerebellar Ataxia or Peripheral Neuropathy .
- Show decline in brain function with accompanying fatigue or Fibromyalgia.
- Have Gluten-Reactivity and/or Dairy-Sensitivity contributing to neurological dysfunction.

 **CYREX**[®]
Laboratories
www.JoinCyrex.com

NEUROLOGICAL AUTOIMMUNE REACTIVITY SCREEN™

